

## Kerraped Boot Order Form – for ordering prescriptions from GP practices

Kerraped boots should only be ordered for those patients who are in compression bandaging which will reduce mobility with a detrimental effect on wound healing and who are at risk of falls from mobilising in unsuitable footwear. Boots are guaranteed by the manufacturer for three months, but may last longer depending on type of use and frequency, and should be measured according to the size guide below to ensure that there is no waste from ordering the wrong size. Patients who are housebound would not normally be expected to require this.


Patient Name	
Patient Date of Birth	
Patient Address	
Name of requesting DN	
Contact number	
Date	
Expected duration of treatment	

Kerraped Size	Shoe Size	Foot width (maximum)	Dressing circumference
Small	UK 2-5.5	7.6 – 8.6cm	23.8 – 27.1cm
Medium	UK 6-7.5	8.8 – 9.2cm	27.6 – 29cm
Large	UK 8-10	9.4 – 10cm	29.5 – 31.4cm
Extra Large	UK 10-13	10.1 – 10.9cm	31.8 – 34.2

How to measure

1. Measure the widest part of the foot without dressings on
2. Measure the circumference of the foot around the widest point with dressings on

How and where to measure:



	Left foot	Right foot
Measurement 1 (cm)		
Measurement 2 (cm)		
Size of boot requested		
Number requested <i>If for both feet please state, are ordered singly and can be used on either foot</i>		